

DigiCAM 2.0

mit USB-LED-Leuchtgriff "Helios 2" (Abb. siehe ACC 21)
with USB-LED-Light handle "Helios 2" (Foto see ACC 21)



Hochwertige Endoskopie-Digitalkamera
mit USB-Kabel für Anschluss an Laptop oder PC,
mit Patientenverwaltung, Bild- und Videospeicher,
geeignet für die Arztpraxis!

Set bestehend aus:

Art.-Nr. :	Bezeichnung :
190.0030	Kamera-Handstück 1/3" CCD-Chip, 440.000 Pixel C-mount-Gewinde
190.0030.K	Kamerakabel USB 2.0 / 2m
190.0030.T	CD mit Treiber für Windows XP/Vista und MedicalEDA 30 Tage Demoversion
190.0030.TK	Mini-Transportkoffer, blau
190.0030.X	Ablagehalter für Handstück
190.0043.PA	TV-Adapter, F = 23 mm
103.3015 + 110.9025	LED-Griff "Helios 2" + USB-Kabel Adapter für Endoskop
190.0031	Software-Lizenz MedicalEDA (für zeitlich unbegrenzte Nutzung)

(Laptop nicht enthalten !)

High-Quality Digital-Camera for Endoscopy
with USB-Cable for connection to Laptop or PC,
with patient administration, image and video capture
suitable for the doctor's office !

Set consisting of :

Ref.-No. :	Description :
190.0030	Camera hand piece 1/3" CCD-Chip, 440.000 Pixel C-mount-Thread
190.0030.K	Camera cable USB 2.0 / 2m
190.0030.T	CD with driver for Windows XP/Vista and MedicalEDA 30 days demo version
190.0030.TK	Mini transport suitcase, blue
190.0030.X	Holder for hand piece
190.0043.PA	TV-Adapter, F = 23 mm
103.3015 + 110.9025	LED-Handle "Helios 2" + USB cable Adapter for Endoscope
190.0031	Software License MedicalEDA (for unlimited use)

(Laptop not included !)