

# FieCAM **CCD-Kamera, 1-Chip, Standard Ausführung** CCD-Camera, 1-Chip, Standard model



	mit/with USB	ohne/without USB
PAL	190.0045.P.USB	190.0045.P
NTSC	190.0045.NT.USB	190.0045.NT

### Digitale 1-Chip-CCD-Farbkamera, bestehend aus:

- Steuereinheit(CCU-Camera Control Unit);
- wasserdichtem Kamerakopf (mit C-mount-Gewinde)
- Kamerakopf-Kabel (NICHT wasserdicht)
- USB-Stick (128 MB) mit integrierter Software; -Fußschalter
- Optional: USB-Funktion, erstellt digitale Bilder über Fußschalter, Speicherung auf USB-Stick

### Technische Daten:

Spannung: 110/115 or 230/240 VAC, 50-60 Hz ; 470 TV-Linien;  
 Shutter: max. 1/2.000.000 sec. ; Videoausgänge 1 x BNC, Y/C (S-VHS); 0,9 lux, f = 1,4,  
 Netzkabel, Bedienungsanleitung, CE, Maße: BxHxT = 215 x 75 x 250 mm, 2,5 kg  
 Optionen: - mit Window-Funktion; - Weiß-Abgleich, - Anti-Moirré-Filter

### Digital 1-Chip-CCD-Colour-Camera, consisting of:

- CCU-Camera Control Unit; -Camera-head (waterproof, with C-mount-thread)
- Camera-head-cable (NOT waterproof);
- USB-Stick (128 MB) with integrated Software; -footswitch
- Optional: USB function, captures digital images over footswitch, saves to USB-Stick

### Technical Data:

Power requirements: 110/115 or 230/240 VAC, 50-60 Hz ; 470 TV-Lines;  
 Shutter: max. 1/2.000.000 sec.; Video output: 1 x BNC, Y/C (S-VHS); Main cable;  
 0,9 lux, f = 1,4; - Service manual, CE, Dim.: mm: 215 x 75 x 250 mm, 2,5 kg  
 Options: - with window function; - white balance; -Anti-Moirré-Filter

