

# LED-Beleuchtungsgriffe für Endoskope


## LED-Illumination handles for Endoscopes

Keine durchgebrannten Glühbirnen mehr !  
Keine unhandlichen Lichtleitkabel mehr !  
Stattdessen absolute Bewegungsfreiheit  
durch innovative Endoskop-Beleuchtungen !

- geeignet für Endoskope bis Ø 6 mm
- Lebensdauer bis 100.000 Betriebsstunden
- passend für Wolf-Ansatz am Endoskop  
(ausser Helios 2 , siehe unten)

No more burned lightsource bulbs !  
No more unhandy light guide cables !  
Instead now absolute mobility through  
your new Endoscope-Illumination:

- applicable for endoscopes up to Ø 6 mm
- Lifetime up to 100.000 operating hours
- fix to Wolf-light-connection at scope  
( except Helios 2 , see below)



### Modell/model

### Art.-Nr./Ref-No.

#### LED 103

Power-LED  
Einfaches Modell  
Länge: 6,5 cm  
Batteriensatz für  
200 Std. Brenndauer

#### 103.1003

Power LED  
Plain Model  
length: 6,5 cm  
battery set for  
200 operating hours



#### LED 203

Power-LED  
Einfaches Modell  
Länge: 14 cm  
Batteriensatz für  
Brenndauer 250 Std.

#### 103.2003

Power LED  
Plain Model  
Length: 14 cm  
battery set for  
250 operating hours



#### HELIOS 2

Hochleistungs-LED  
Exquisites Modell  
Länge: 7 cm, Kabel 2 m  
Akkuladung für 4 Std. Betrieb  
Farbtemperatur 5.500° Kelvin  
Adapter für Endoskope von:

#### 103.3010

High-Output-LED  
Exquisite model  
Length: 7 cm, cable 2 m  
Loaded accu for 4 operat. hours  
Color temperature 5.500° Kelvin  
Adapter for endoscopes of:

|                          |          |
|--------------------------|----------|
| R. Wolf                  | 110.9022 |
| Olympus                  | 110.9025 |
| ACMI<br>British Style    | 110.9024 |
| K. Storz<br>Fiegert-End. | 110.9025 |

|                          |          |
|--------------------------|----------|
| R. Wolf                  | 110.9022 |
| Olympus                  | 110.9025 |
| ACMI<br>British Style    | 110.9024 |
| K. Storz<br>Fiegert-End. | 110.9025 |

