




Automatik-Trokarhülse Automatic trocar sleeve



Automatik-Trokarhülse, mit Insufflationshahn und manueller Klappenöffnung

Automatic trocar sleeve, with insufflation stop cock and manual opening of trap door valve

Ø in mm	Nutz- länge Length	Trokar- hülse Sleeve	Trokare / trocars	
			Dreikant Pyramidal	Konisch Conical
				
5,0	92 mm	213.0050	213.0050.D	213.0050.K
5,5	95 mm	213.0055	213.0055.D	213.0055.K
10,0	97 mm	213.0010	213.0010.D	213.0010.K
11,0	101 mm	213.0011	213.0011.D	213.0011.K
12,5	98 mm	213.0012.5	213.0012.5D	213.0012.5K

Zubehör für Trokarhülsen Accessories for trocar sleeves



für Trokarhülsen mit
for trocar sleeves with

5 + 5,5mm Ø	210.0005
10mm Ø	210.0010
11mm Ø	210.0011
12,5mm Ø	210.0012

Fixier- Spiralhülse zum perkutanen Eindrehen und Fixierung von Trokarhülsen.
Fixation sleeve spiral type for percutaneous screwing and fixation of trocar sleeves