

## Sicherheitstrokar Safety trocar



Sicherheitstrokar wiederverwendbar, schneidend,  
mit automatischem Schneidenschutz, zerlegbar  
passend zu Standard Trokarhülsen

Safety-trocar reuseable, cutting, with  
automatic edge cover, dismountable,  
fitting to standard trocar sleeves



Ø	Ref.-No.
5,0 mm	210.0504
5,5 mm	210.5504
10 mm	210.1004
11 mm	210.1105
12,5mm	210.1204
15 mm	210.1504

## Reduzierer - Reducer

Reduzier-Hülse  
Reducer sleeve



	Ref.-No.		Ref.-No.
5 mm - 4 mm	210.0506	15 mm - 5,0/5,5 mm	210.1509
10 mm - 5,0/5,5 mm	210.1008	15 mm - 10 mm	210.1508
11 mm - 5,0/5,5 mm	210.1108	20 mm - 5,5 mm	210.2055
12,5 mm - 5,0/5,5 mm	210.1208	20 mm - 10 mm	210.2010
12,5 mm - 10 mm	210.1209	20 mm - 15 mm	210.2015

## Dilatation Dilation

Dilatations-Satz mit Führungsstab  
Dilation set with guide probe



5mm/10mm	210.1006
5mm/11mm	210.1107
10mm/15mm	210.1506
11mm/15mm	210.2006
11mm/20mm	210.2009