














# Modulare bipolare Koagulationsfaßzange

## Modular bipolar coagulation forceps



**Ø 5 mm, Nutzlänge / Working length: 360 mm**

	<b>Tubensterilisation</b>	Tube sterilization	<b>230.5360</b>
	<b>Faßzange, gefenstert</b>	Grasping forceps	<b>230.5361</b>
	<b>Maul: 3 mm breit, glatt</b>	Jaw: 3 mm wide, smooth	<b>230.5362</b>
	<b>Maul: 3 mm breit, gerieft</b>	Jaw: 3 mm wide, grooved	<b>230.5363</b>
	<b>Maul: 1,5 mm breit, gerieft</b>	Jaw: 1,5 mm wide, smooth	<b>230.5364</b>
	<b>Maul: 1,5 mm breit, glatt,</b>	Jaw: 1,5 mm wide, grooved	<b>230.5365</b>
	<b>Mikro-Fasszange</b>	Micro grasping forceps	<b>230.5366</b>
	<b>Ersatzteile und Zubehör</b>	Spare parts and accessories	:
	<b>U-Griff</b>	U-Handle	<b>230.5360.G</b>
	<b>Schaftrohr</b>	Sheath tube	<b>230.5360.S</b>
	<b>Ringgriff</b>	Ring-Handle	<b>230.5370.G</b>
	<b>Innenteile: Art.-Nr. wie oben + ".I"</b> (Beispiel 230.5360.I , 230.5361.I , ...)	Inner parts: Ref.-No. as a.m. + ".I" (f.ex. 230.5360.I , 230.5361.I , ...)	
	<b>Bipolar-Kabel</b> <b>Länge 3 m, zu Geräten von:</b>	<b>Bipolar cables</b> <b>length 3 m, to devices of</b>	:
	<b>Storz, Wolf, Erbe, Wisap</b>	dto.	<b>230.5350</b>
	<b>Berchtold, Martin</b>	dto.	<b>230.5350.M</b>
	<b>USA-Devices</b>	dto.	<b>230.5350.A</b>