


Tubensterilisation, ein Einstich

Tube sterilization, single puncture

OP-Laparoskop Ø 10 mm
Blickrichtung 0°
Arbeitskanal Ø 6 mm
Nutzlänge 268 mm
autoklavierbar, inklusive
abnehmbarem Einführungshahn

Parallel-Einblick
Parallel eye piece

201.1000.A



Operating laparoscope Ø 10 mm
direction of view 0°
working channel Ø 6 mm
working length 268 mm
autoclavable, inclusive
detachable instrument port

45° - Schräg-Einblick
Angled eye piece

201.1100.A



Einführungshahn / Instrument port:

-mit Standard-Ventil
-with standard-valve

200.1101



-mit Automatik-Ventil
-with automatic valve

201.5501



Ring-Applikator, ohne Ladevorrichtung
Ring-Applicator, without loading device

Ø 6 mm

250.0600

Ladevorrichtung, autoklavierbar (Konus + Schieber)
Loading device, autoclavable (Cone + Slide)

(250.0601.C / 250.0601.S)

250.0601