


## Naso-Pharyngo-Laryngo-Fiberscopes " SILVER-LINE " + " SILVER-LINE SLIM "



Ø 3,4 mm " SILVER-LINE " 800.3800  
Ø 2,8 mm " SILVER-LINE SLIM " 800.2801



800.3800.K

### Technische Daten:

- Nutzlänge 30 cm
- Abwinkelung auf-/abwärts 150°
- Blickrichtung 0°
- Bildwinkel 70°
- Schärfebereich 1-50 mm
- Gas sterilisierbar
- einlegbar

### - inklusive Standard-Zubehör:

- Transportkoffer (1)
- ETO-Kappe für Druckausgleich (2)
- Dichtigkeitstester (3)
- Gebrauchsanweisung

### Technical Data:

- Working length 30 cm
- Angulation up-/downward 150°
- Direction of view 0°
- Field of view 70°
- Depth of field 1-50 mm
- Gas sterilisation
- soakable

### - inclusive Standard-Accessories:

- Transport suitcase (1)
- ETO-Cap for pressure compensation (2)
- Leakage tester (3)
- Instructions manual



800.5009



800.5010

**Steriler Einweg-Überzug für 800.2801 - VE = 50 Stück 800.2801.X**  
Sterile single use cover for 800.2801 - PU = 50 pieces

**Steriler Einweg-Überzug für 800.3800 - VE = 50 Stück 800.3800.X**  
Sterile single use cover for 800.3800 - PU = 50 pieces

